# **CEL**signia

Hyperactive c-Met and ErbB signaling detected in a sub-group of ovarian patient tumors: Patient sub-group may benefit from c-Met and pan-HER combination therapy

Stefano Rossetti, Catherine Kuzmicki, Aaron Broege, Joanna Sabat, Salmaan Khan, Ian MacNeil, Brian Sullivan, Lance Laing | Celcuity, 16305 36th Ave N, Suite 100, Minneapolis, MN 55446, USA

(A) Table showing HER receptors expression

(A) Table showing HER receptors expression (FACS stain index) and CELsignia responses: HER1 signaling (EGF activation), HER3 signaling (NRG activation), HER2 signaling (CELsignia)

HER2 Test score), c-Met (HGF activation), The

breast cancer cells lines MDA-MB-231 (with HER2 expression "0" by IHC), MDA-175 (with

IER2 "+1" by IHC) and SKBR3 (with HER2 "+3 by IHC) are shown as comparison

# Background

Sub-groups of breast, lung, and colon cancer patients with genetic variants associated with dysregulated ErbB (HER) signaling benefit from treatment with ErbB targeted therapies. No sub-groups of ovarian cancer (OxC) patients with actionable ErbB genetic variants have yet been found. Measurement of biological factors other than genetic variants, such as ErbB signaling activity, may identify ovarian cancer patients likely to benefit from ErbB targeted therapies. We have previously reported studies detecting dysregulated ErbB and c-Met signaling activity in 20%-25% of HER2-negative breast cancer patient tumors using the CEJsina Multi-Pathway (MP) Signaling Function Test. To determine whether ErbB or c-Met signaling dysregulation is involved in ovarian cancer, the CELsignia MP test was adapted to analyze ovarian tumor cells

1) characterize c-Met and ErbB family signaling activity in ovarian patient tumor cells and ovarian tumor cell lines; 2) evaluate in vivo response to pan-HER and c-Met inhibitors in an ovarian tumor xenograft model

### Methods

Cell Lines: The ovarian cancer cell lines used in this study (listed in Table 1) were maintained according to ATCC recommendations and authenticated by ATCC.

Tissue Specimens: A training set of de-identified ovarian cancer tissue specimens was obtained from 19 patients patient characteristics are listed in Table 2). Methods for tissue extraction and primary cell culture were based on (Laing et al. 2016, 2017).

	OvC type	Sample	Age	Ethnicity	Primary	Removal	Туре	Grade
CaOV-3	High Grade Serous	C1113.01.38	53	Caucasian	Ovary	Ovary	Serous	High grad
CaOV-4	High Grade Serous	C1468.01.18	65	Caucasian	Ovary	Ovary	Serous	High grad
COV362	Endometriaid/HGS	C1497.01.10	58	Caucasian	Ovary	Ovary	Serous	High grad
HEY-T30	High Grade Serous	C1504.01.10	58	African Amer.	Ovary	Ovary	Serous	High grad
NURAMOCHI	High Grade Serous	C1508.01.13	45	Caucasian	Ovary	Ovary	Endometriold	Grade 2
MIS-OV	Serous Cystadenocarcinoma	C1620.01.21	42	Caucasian	Ovary	Omentum	Serous	High grad
OV-90	Adenocarcinoma	C1774.01.11	45	Caucasian	Ovary	Omentum	Serous	High grad
OVCAR-3	High Grade Serous	C1928.01.4	67	African Amer.	Ovary	Ovary	Serous	High grad
OVCAR-4	High Grade Serous	C1944.01.3	52	Caucasian	Ovary	Ovary	Serous	High grad
OVKATE	High Grade Serous	C1966.01.47	58	Caucasian	Ovary	Ovary	Serous	High grad
OVMANA	Clear Cell Adenocarcinoma	C2032.01.4	78	Caucasian	Ovary	Omentum	Serous	High grad
OVSAHD	High Grade Serous	C2040.01.9	49	Asian	Ovary	Ovary	Serous	High grad
SK-OV-3*	Serous Cystadenocarcinoma	C2060.01.1	59	Caucasian	Ovary	Ovary	Serous	High grad
TOV-112D	Endometriald	C2064.01.1	70	Caucasian	Ovary	Omentum	Serous	High grad
TOV-21G	Clear Cell Adenocarcinoma	C2095.01.1	63	Caucasian	Ovary	Ovary	Serous	High grad
UACC-1598	Serous Cystadenocarcinoma	C2332.01.10	55	Caucasian	Ovary	Omentum	Serous	High grad
UACC2727	Serous Carcinoma	C2408.01	59	Caucasian	Ovary	Ovary	Serous	High grad
UW01.289	Ovarian carcinoma	C2427.01.5	48	Caucasian	Ovary	Omentum	Serous	High grad
Reported to have Her2 amplification		C2457.01.1	64	Caucasian	Quary	Overy	Secola	High grad

Flow cytometry: Flow cytometry on disaggregated tissue and cultured cells was performed on Agilent Novocyte 3005. Extractular staining with indicated antibodies was performed by standard procedures. DNA staining for DNA index calculation was performed with FxCycle violet.

Sequencing: Genomic DNA was extracted from tissue or cells and provided to contract research organizations for targeted next-generation sequencing to an average depth of 1000x. Sequencing was targeted to panels of 500 – 1000 genes frequently mutated in solid turnors with enrichment (Nimblegen SecQar or Novogene NovOVPM 1.0). Reads were aligned to the reference genome with BWA. Tumor vs. normal comparisons were performed with VarScan or Mutect<sup>™</sup> programs

Xenografts: 40 female NSG™ mice were injected with two million OVCAR-4 cells. Mice were randomly assig to either a control group that received Captiol excipient or a treatment group that received either neratini teportinib, on reactinib and tepotionib for 16 days.

CELsignia analysis: Real-time live cell response to specific ErbB and c-Met agonists (3 nM NRG1b, 0.3 nM EGF, and 4 ng/mL HGF) with or without an antagonist (2C4, a HER2 dimerization inhibitor; neratinib, a pan-HER inhibitor; or tenotinib, a c-Met kinase inhibitor) was measured and quantified using an xCELLigence impedance biosensor or reporting, a criver kilose financio y las mesoareo an uparather using an Accurgence indexessor (Aglent Technologies). Signalia estivity above an unevolutive stabilistica cui of value (25) signalia purify vas used to identify abnormal levels of EGFR, HER2, and c-Met signaling activity. A HER2 Signaling Score, reflecting NRG and EGF responses that can be attenuated by 2C4, was acciulated as was previously described (Laing et al. 2016, 2017). Figure 1. CELsignia Analysis Uses Biosensor to Quantify Signaling Activity in Real-Time in Live Cells • CFI signia analysis leverages connections between cell adhesion impedance, and cell signal transduction Live cells are attached via ECM to a microelectrode on the bottom of a 96-well impedance biosensor plate. Additionally, the tumor epithelial cells in the wells form adhesion-based gap junctions

The cells attached to the biosensor impede the flow of electrons when mV AC current is applied and changes in impedance (mO) are recorded





Results

# Figure 2. CELsionia analysis identifies OvC cells lines with hyperactive c-Met and ErbB signaling



Δ



These results show that : 1. 3 of 18 OvC cell lines (OVMANA, OVACR3, OVCAR4) have hyperactive c-Met and/or ErbB signaling 2. Hyperactive ErbB signaling does not correlate with HER1-4 expression 3. Neratinib significantly attenuates NRG and EGF response 4. Tepotinib significantly attenuates HGF signaling respons



These results show that neratinib and tepotinib induce significant tumor reduction in an OVCAR4 xenograft model found by the CELsignia test to have hyperactive ErbB and c-Met signaling.

Figure 3: Neratinib and tepotinib reduce xenograft tumor growth in OvC cells with hyperactive c-Met

### Figure 4: Establishment of ovarian cancer primary cultures for CELsionia analysis



Digest Culture

1.29 1.38 1.70 1.09

1.04 1.15 1.66

1.56 2.00 1.59 1.03

0.93

1.68 0.98 1.21 1.75

1.88

1.54 1.64 1.72 1.84

1.37 1.51 1.31 1.78

7.9 4.7

13 11 42 23

cells (EpCAM-/Podoplanin+) and Digest Culture Digest Culture C1113.01 C1468.01 C1497.01 21504.01 21508.01 (C) Genetic analysis confirms that cultured cells are enriched for NA 8.9 4.0 4.5 mutations present in cancer tissue 50.1% (selected mutations are shown). 7.0 9.8 (D) Screening of 19 OvC natient (D) Screening of 19 OVC patient specimens (see Table 2) showing that primary cultures are enriched for epithelial (EpCAM+) cells with a DNA index similar to what found in the 1928.01 digest tissue. FACS reference ranges for HER2 expression: SKBR3 (IHC 3+ 408); MDA-175 (IHC 1+; 66); MDA-231

## These results show that

(IHC 0; 13).

(B) Representative FACS analysis of an initial disaggregated tissue and the resulting cultured cells showing similar

% of epithelial (EpCAM+) and stromal

similar hyperploid DNA content

1. Primary cultures are enriched for epithelial cancer cells present in the original OvC tissue 2. HER2 is expressed at low levels in all specimens tested.

#### Figure 5: CELsignia analysis identifies OvC patients with hyperactive c-Met and ErbB signaling

celcuity

www.celcuity.com

500



(A) Table showing HER recentors expression (FACS stain index) and CEI signia responses: HER1 signaling (FGE activation) (4) Table showing mark receipus expression (rAcc 3 stain marks) and CLSIgnin responses. Track agring (CEF HER3 signaling (NRG activation), HER2 signaling (CELSignal HER2 Test score), c-Met (HCF activation). FACS refi ranges for Her2 expression: SKBR3 (HC 3+; 408); MDA-175 (HC 1+; 66); MDA-231 (HC 0; 13). (B) Plot showing lack of correlation between HSF score and HER2 expression (FACS stain index

These results show that : 1. 2 of 19 OvC patient samples have hyperactive c-Met and ErbB signaling. 2. Hyperactive ErbB signaling does not correlate with HER1-4 expression.

# Summary of Results

The CELsignia Multi-Pathway Signaling Function Test identified 3 of 18 (17%) OvC cell lines and 2 of 19 (11%) OvC patients with hyperactive c-Met and/or ErbB signaling.

 FACS and sequencing analyses confirmed that OvC primary cultures maintained the epithelial c phenotype present in the original tissue. Hyperactive ErbB signaling did not correlate with HER1-4 expression

Neratinib and tepotinib significantly attenuated ErbB signaling and c-Met signaling in the CELsignia

In vivo studies showed that mouse xenograft tumors developed from an OvC cell line with hyperactive ErbB and c-Met signaling significantly decreased in size after treatment with a combination of tepotinil

# Conclusions

group of ovarian cancer patients have abnormal ErbB a treatment with a combination of ErbB and c-Met inhib e of this patient sub-set to combined c-Met and pan-H

## References

 Laing, et al. A functional signal profiling test for identifying a subset of HER2-negative breast cancers with abnormally amplified HER2 signaling activity. Oncotarget. 2016. 7(48):78577-78590. Laing, et al. Development of a test that measures real-time HER2 signaling function in live breast cancer cell lines and primary cells. BMC Cancer 17:199.

This presentation is the intellectual property of Celcuity Inc. Contact Brian Sullivan at BSullivan@Celcuity.com for permission to reprint and/or distribute.